

Dealer _____

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 May 2024

Key ☐ standard ☐ options


Dampened backrest uses gas springs to compensate for rapid involuntary body movements.

DAMPENED backrest ALTHEA.M

Alum. foldable frame

- L** Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46 ☐48 ☐50
☐52 ☐56 ☐60 (including cover for footplate according to setup)
- P** Seat depth ☐36 ☐38 ☐40 ☐42 ☐44 ☐46
 Adjustable up to 6cm (max.50 cm) ☐48 ☐50
- HS** Back-rest height ☐29 ☐31 ☐33 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45
☐47 ☐49 ☐51 ☐55 ☐60
- HA** Front seat/floor height ☐42 ☐43 ☐44 ☐46 ☐48 ☐50 ☐52
☐54
- HP** Rear seat/floor height ☐38 ☐40 ☐42 ☐44 ☐46 ☐48 ☐50
- DP** Foot-rest/seat distance ☐34 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45 ☐47 ☐49 ☐51
 (from 17 to 32 cm foot-rest directly assembled on the frame)

cm



- Frame colour** ☐ Matte black
☐ Metallic blue ☐ Pearl red ☐ Metallic grey
☐ Orange ☐ Emerald green ☐ Deep sky light blue
☐ (see "Colours" catalogue) Neon colour

①

Upholstery

- ☐ Tension adj. nylon backrest upholstery
☐ Tension adj. BREATH.backr.uph.
☐ Without any backrest upholstery
☐ With divided seat (joint by velcro, def.)
☐ With one entire seat (not depth adjustable)

②


Upholstery colour

(Coloured uphol. only with Tension. adj. nylon)


- ☐ Black
☐ Red * ☐ Orange *
☐ Light blue * ☐ Green *
☐ Cornflower blue *

③


- ☐ Not Detachable foot-rest
☐ Dampened not det. foot-rest
 (needs single dampened footplate with fastening block)




- ☐ Detachable foot-rest 82° angle
 (DP 41 cm)
 (82° footrests do not allow any 175 front castors)




- ☐ Detachable foot-rest 68° angle
 (DP 41 cm)




- ☐ Dampened det. foot-rest 68° angle
 (needs single dampened footplate with fastening block)



- ☐ Elevating foot-rest
 (detachable)



- ☐ Detachable Desk
 armrest h 22 cm



- ☐ 27 cm high from the seat


Elbow-rest

- ☐ Classic (30 cm long)
☐ Straight (26 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

⑧

- ☐ Plastic separate





- ☐ Cover f/footplates *
☐ Foot-holding belts (W50max)
☐ Heel stops (W50max)
 (* eventually included for W 52 to 60)

⑤

Single aluminium profile tilt adj.footplate:

(not available with elevating footrest)



- ☐ forward ☐ backward



- ☐ Foot-holding belts

Separate aluminium profile tilt adjustable footplates:

- ☐ forward ☐ backward




- ☐ Foot-holding belts (W50max)

⑥

Single DAMPENED footplate aluminium profile, w/fastening block (needed in case of dampened footrest)


- ☐ forward
☐ backward



- ☐ Foot-holding belts

⑦

- ☐ Detachable "U"
 armrest h 22 cm



- ☐ 27 cm high from the seat


Elbow-rest

- ☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX


note: see also Hemi modular armrest ord.form

⑨

Height adjustable detachable armrest (only closed)



- ☐ Desk model



- ☐ "U" model

Elbow-rest

- ☐ Classic (30 cm long)
☐ Classic (36 cm long)
☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

⑩

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 May 2024

DAMPENED backrest
ALTHEA.M

Dealer _____

Rear wheels ⑪

- ☐ 24" 36 spokes ☐ 12" solid*
☐ 22" 36 spokes ☐ 12" pneumatic*
☐ 20" 36 spokes (* N.A. if HP=42 or HP≥46cm if HP44: only curved rear frame)
☐ 24" MTB citybike tyre
☐ 24" MTB offroad tyre

- ☐ Pneum. Plugged (24x1.3/8") **Tyres**
☐ Plugged ☐ Smooth Pneum. high press.
☐ Solid grey 24x1.3/8"
☐ Solid KIK 24" black
☐ 24" Pneum. Marathon Plus

Pushrims ⑫

- ☐ Aluminium silver
☐ Chromium-plated steel 24"
☐ Aluminium dark
☐ Rubber coated
☐ Carbolife CURVE L 24"(*)
☐ OVAAL silver alum. (*)
☐ OVAALS dark alum. (*)
☐ OVAALG dark, w/silic. insert (*)
☐ Carbolife GEKKO 24" (*)
☐ Silicon covers 24"
 (*) not available with MTB wheels

Spokeguards ⑬



- ☐ Transparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)

Front castors ⑭

- ☐ 100 (4"x1.1/4) solid - plastic rim
☐ 100 (4"x1.1/4) solid - alum. rim
☐ 125 (5"x1.1/4) solid - plastic rim
☐ 125 (5"x1.1/4) solid - alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid
☐ 175 (7") solid **
☐ 175 (7") pneumatic **
 (**) only if HA ≥48cm and with 68°detach.footrests

STRAIGHT Rear Frame

- ☐ With small wheels for narrow passages (default) ☐ Without any small wheels for narrow passages




Unbalancing system

- ☐ Right
☐ Left
 (Only with straight rear frame)

Anti-tip device ⑮

- ☐ Right
☐ Left
 (Only with straight rear frame)

CURVED Rear Frame Unbalancing system



- ☐ Right
☐ Left
 (Only with curved rear frame)

Revolving anti-tip device


- ☐ Right
☐ Left
 (Only with curved rear frame)

(With curved rear frame it is not possible to assemble any small wheels for narrow passages)

Single drive ⑯

(only with 24" or 22" rear wheels)

- ☐ Right
☐ Left




Push handles ⑰

- ☐ Standard ☐ Height adjustable
 (max L = 48 cm)
 (Compulsory with HP44cm or lower)




Extended wheel plate ⑱



(The rear wheels are assembled with COG -2cm)

Lateral Supports, swing away ⑲




- ☐ Right
☐ Left

Specify the height "h"
 from center of Swing Away to seat

Head-rest ⑳

Padded




- ☐ Fixed
☐ Height adjustable

Orbital


- ☐ Anatomic



- ☐ Moulded



- ☐ With lateral supports




Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm

h = cm d = cm

Table ㉑

- ☐ Transparent
☐ RxOLx Hemiplegic transp.
☐ Plastic grey
☐ Soft padded

Specify the fastening:
☐ single ☐ double



(Only with arm-rest) if L>50cm

Lined cushion in foam polyurethane ㉒



- ☐ 3 cm high
☐ 5 cm high


Brakes ㉓

- ☐ Classic - short handle
☐ Classic - long handle
☐ Drum
☐ Right single brake
☐ Left single brake
☐ Drum+classic

Options

- 25 ☐ "M" velcro strip (hook) on seat
 26 ☐ "M" velcro for VICAIR
 27 ☐ "F" velcro strip (loop) on seat
 28 ☐ Seat/back-rest connection flap
 29 ☐ Chest harness
 30 ☐ Chest belt
 31 ☐ Abductor belt
 32 ☐ Front protective covers
 33 ☒ Stretched bar
 34 ☐ Catheter holder pocket
 (note: only with entire seat)

Abductor ㉔



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm

☐ cm
 Choose or specify the height from seat upholstery

- 35 ☐ R ☐ L Crutch hold./single
 36 ☐ R ☐ L Crutch hold./double
 37 ☐ Pump for pneumatics
 38 ☐ Suitcase holder
 39 ☐ 45° pelvic belt
 40 ☐ 45° padded pelvic belt
 41 ☐ Black backpack
 42 ☐ IV bar + poles
 43 ☐ IV poles

- 44 ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)
 45 ☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
 46 ☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes