

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 May 2024

Dealer


Key

☐ standard

☐ option

**MINISTAR &
 MINISTAR PREMIUM**

L Seat width	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="radio"/> 46 <input type="radio"/> 48 <input type="radio"/> 50 (compulsory stretched bar to be added) (included double cross bar) <input type="radio"/> 52 <input type="radio"/> 54 (compulsory stretched bar to be added) (included double cross bar) <input type="radio"/> 56 <input type="radio"/> 58 <input type="radio"/> 60 (compulsory stretched bar to be added) (included double cross bar)
P Seat depth	<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="radio"/> 44 <input type="radio"/> 46 <input type="radio"/> 48* <input type="radio"/> 50*(*compulsory double crossbar to be added)
HS Back-rest height	<input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 <input type="radio"/> 43 <input type="radio"/> 45 <input type="radio"/> 47 <input type="radio"/> 49 <input type="radio"/> 51 (included stretched bar)
HA Front seat/floor height	<input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="checkbox"/> 52 <input type="radio"/> 54
HP Rear seat/floor height	<input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50
DP Foot-rest/seat distance	<input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 <input type="checkbox"/> 43 <input type="checkbox"/> 45 <input type="checkbox"/> 47 <input type="checkbox"/> 49 <input type="checkbox"/> 51 cm



Frame colour

- ☐ Metallic blue ☐ Pearl red ☐ Metallic grey
☐ Glossy black ☐ Matte black
☐(see "Colours" catalogue)
☐ Neon green ☐ Neon orange

①

Upholstery

- ☐ Nylon (L50 max)
☐ Tension adj. nylon back-rest
☐ Tension adj. breathable back-r.
☐ Without any back-rest upholstery

②


Upholstery colour

(nylon only)


- ☐ Black
☐ Red ☐ Orange
☐ Light blue ☐ Green
☐ Cornflower blue

③

☐ **Not-detachable foot-rests**




☐ **Detachable foot-rests**



☐ **Elevating foot-rests (detachable)**

④


☐ **Plastic separate**
 (only L = from 38 to 46 cm)



- ☐ Foot-holding belts
☐ Heel stops

⑤

Separate aluminium profile tilt adjustable foot-plates:




- ☐ Foot-holding belts


⑥

Single aluminium profile tilt adjustable foot-plate:
 (not available with elevating foot-rest)

☐ forward



☐ backward




- ☐ Foot-holding belts

⑦

Side guards
 (not available w/reclining backrest)

☐ Plastic



☐ Aluminium
 (same colour)



☐ Carbon fibre




☐ Plastic round cover



⑫

☐ **Detachable Desk arm-rest - h 22 cm**



- ☐ Open* with side guards (only with HP > 42 cm)
☐ Closed
☐ 27 cm high from the seat
☐ Tip-up*
 *(not available with reclining backrest)


Elbow-rest

- ☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

⑧

☐ **Detachable Sport arm-rest - h 22 cm***
 *(not available with reclining backrest)



- ☐ Open with side guards (only with HP > 42 cm)
☐ Closed
☐ 27 cm high from the seat
☐ Tip-up


Elbow-rest

- ☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod. armrest ord.form

⑨

☐ **Detachable "U" arm-rest - h 22 cm**



- ☐ Open* with side guards (only with HP > 42 cm)
☐ Closed
☐ 27 cm high from the seat
☐ Tip-up*
 *(not available with reclining backrest)

Elbow-rest

- ☐ Straight (26 cm long)
☐ Classic (36 cm long)
☐ Tubular
☐ Hemiplegic RX
☐ Hemiplegic LX


note: see also Hemi mod. armrest ord.form

⑩

Height adjustable arm-rest (only closed)



- ☐ Desk model
☐ Tip-up*
 *(not available with reclining backrest)



☐ "U" model

- ☐ Tip-up*
 *(not available with reclining backrest)

Elbow-rest

- ☐ Classic (30 cm long)
☐ Classic (36 cm long)
☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod. armrest ord.form

⑪

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 May 2024

MINISTAR &
MINISTAR PREMIUM

Dealer _____

Rear wheels

12

- ☐ 24" (36 spokes) ☐ Light 24"
☐ 22" (36 spokes)
☐ 20" (36 spokes)
☐ 24" Mountain-bike w/Citybike tyre
☐ 24" Mountain-bike w/Offroad tyre

Tyres

- ☐ Pneum. Plugged (x 1"3/8)
☐ Solid grey x 1"3/8
☐ Pneum. high press. Plugged
☐ Pneum. high press. Smooth
☐ 24" Pneum. Marathon Plus

Pushrims

13

- ☐ Aluminium silver
☐ Chromium-plated steel 24"
☐ Aluminium dark
☐ Rubber coated
☐ Carbolife CURVE L 24"(*)
☐ Carbolife GEKKO 24" (*)
☐ Silicon covers 24"

(*) not available with MTB wheels

Spoke guards

14

- ☐ Transparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)

Front castors


15

- ☐ 125 (5") solid
☐ 150 (6") solid
☐ 175 (7") pneumatic
☐ 175 (7") solid


Rear frame

16


- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)



- ☐ STRAIGHT rear frame without any small wheels for narrow passages




- ☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels



(Only if the rear seat/floor height HP is between 44 and 48 cm.
 Not selectable in case of reclining backrest)

- ☐ CURVED rear frame




(It is not possible to assemble the small wheels for narrow passages)

Unbalancing system

17




(Only with straight rear frame)



(Only with curved rear frame)

Anti-tip device

(Only with straight rear frame)



- ☐ Right
☐ Left

Revolving anti-tip device

(Only with curved rear frame)




- ☐ Right
☐ Left

Single drive


18

- ☐ Right
☐ Left



Extended wheel plate

19



Compulsory with reclining back-rest

Pushing handles


20

- ☐ Standard



- ☐ Hinged back-rest

(no stretched bar)
 (max L = 48 cm)
 (Not with reclining back-rest)
 (need tension adj. backrest)




- ☐ Height adjustable
 (max L = 48 cm)
 (need tension adj. backrest)



Reclining back-rest

21

- ☐ from 90° to 150° with fixed pushing handles
☐ from 90° to 150° with height adjustable pushing handles (included)
 (max L = 44 cm)




Compulsory:

- closed arm-rest
 - stretched bar
 - extended wheel plate
 - tension adj. backrest upholstery

Abductor

22




- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm

- ☐ cm
 Choose or specify the height from seat cover

Lateral supports

23




- ☐ Right
☐ Left

Specify the height "h"
 (need tension adj. backrest)

Head-rest

24

Padded




- ☐ Fixed
☐ Height adjustable

Orbital

- ☐ Anatomic



- ☐ Moulded



- ☐ With lateral supports



Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = cm d = cm

Table

25

- ☐ Transparent
☐ Rx ☐ Lx Hemiplegic transp.
☐ Plastic grey
☐ Soft padded

Specify the fastening:

- ☐ single ☐ double



(Only with arm-rest) if L > 50cm

Lined cushion

28

in foam polyurethane



- ☐ 3 cm high
☐ 5 cm high

Brakes

26

- ☐ Classic - short handle
☐ Classic - long handle
☐ Right single brake
☐ Left single brake
☐ Drum
☐ Drum + classic

Options

- 29 ☐ "M" velcro strip (hook) on seat
 30 ☐ "F" velcro strip (loop) on seat
 31 ☐ Seat/back-rest connection flap
 32 ☐ Chest harness
 33 ☐ Chest belt
 34 ☐ Abductor belt
 35 ☐ Front protective covers
 36 ☐ Stretched bar
 37 ☐ Catheter holder pocket
 38 ☐ Double crossbar

€ 115,00

Anodized set

27

(Rims & hub of LIGHT rear wheels and forks)

- ☐ Black
☐ Silver
☐ Orange
☐ Light blue
☐ Red
☐ Cyclamen purple

Options

- 39 ☐ R ☐ L Crutch hold./single
 40 ☐ R ☐ L Crutch hold./double
 41 ☐ Pump for pneumatics
 42 ☐ Suitcase holder
 43 ☐ 45° pelvic belt
 44 ☐ 45° padded pelvic belt
 45 ☐ Black backpack
 46 ☐ IV bar + poles
 47 ☐ IV poles
 48 ☐ Padded rigid cushion (L52max)
 49 ☐ Padded rigid cushion (L56-60)