

Order N° _____ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____


Key
 standard
 option

PAGE 1/2
 May 2024

Dealer _____

- L** Seat width 34 36 38 40 42 44 46 48 50
52* 56* 60 (including cover for footplate according to setup)
- P** Seat depth 36 38 40 42 44 46
 Adjustable up to 6 cm (max.50cm) 48 50
- HS** Back-rest height 29 31 33 35 37 39 41 43 45
47 49 51 55 60
- HA** Front seat/floor height 38_(fc80max) 39_(fc100max) 40_(fc125max) 42_(fc150max) 45_(fc175max)
 (with MAX size of front castors)
 (NOTE: 175 front castors do not allow 82° detachable footrests)
44_(fc150max) 46_(fc175max)
- HP** Rear seat/floor height 34 36 38 40 42 44 45 46 48
 NOTE: 24" rear wheels (if HP min.38/max 48) 22" (if HP min.36/max 45) 20" (if HP min.34/max 43)
- DP** Foot-rest/seat distance 34 35 37 39 41 43 45
 (from 17 to 32 cm footrest directly assembled on the frame)

ALTHEA19.M
LOWERED set up
& DAMPENED backrest
 alum.foldable frame



cm Dampened backrest uses gas springs to compensate for rapid involuntary body movements.

- Frame colour**
- Matte black
 Metallic blue Pearl red Metallic grey
 Orange Emerald green
 (see "Colours" catalogue)
 Neon orange Neon green

- 1 Upholstery**
- Tension adj. nylon backrest upholstery
 Tension adj. breath.backr.
 Without any back-rest upholstery
 With divided seat (joint by velcro, def.)
 With one entire seat (not depth adjustable)

- 2 Upholstery colour**
- (Coloured uphol. only with Tension. adj. nylon)
- Black
 Red* Orange*
 Light blue* Green*
 Cornflower blue*

- Not Detachable foot-rest**
 Dampened not det. foot-rest
 (needs single dampened footplate with fastening block)
- Detachable foot-rest 82° angle**
 (DP 41 cm)
 (82° footrests do not allow any 175 front castors)
- Detachable foot-rest 68° angle**
 (DP 41 cm)
- Dampened det. foot-rest 68° angle**
 (needs single dampened footplate with fastening block)
- Elevating foot-rest (detachable)**

- Plastic separate**
-
- Cover f/footplates*
 Foot-holding belts (L50max)
 Heel stops (L50max)
 (* eventually included f/L 52 to 60)

- 3 Separate aluminium profile tilt adjustable footplates:**
- forward backward
-
- Foot-holding belts (L50max)

- 4 Single aluminium profile tilt adj.footplate:**
 (not available with elevating footrest)
- forward backward
-
- Foot-holding belts

- 5 Single DAMPENED footplate**
 aluminium profile, w/fastening block
 (needed in case of dampened footrest)
- forward backward
-
- Foot-holding belts

- Detachable Desk armrest h 22 cm**
-
- 27 cm high from the seat
- Elbow-rest**
- Classic (30 cm long)
 Straight (26 cm long)
 Tubular
 Hemiplegic RX
 Hemiplegic LX

- Detachable "U" armrest h 22 cm**
-
- 27 cm high from the seat
- Elbow-rest**
- Straight (26 cm long)
 Straight (36 cm long)
 Tubular
 Hemiplegic RX
 Hemiplegic LX

- 6 Height adjustable detachable armrest**
 (only closed)
- Desk model
-
- "U" model
-
- Elbow-rest**
- Classic (30 cm long)
 Classic (36 cm long)
 Straight (26 cm long)
 Straight (36 cm long)
 Hemiplegic RX
 Hemiplegic LX

Dealer _____

Rear wheels (13)

- 24" 36 spokes 12" solid*
- 22" 36 spokes 12" pneumatic*
(N.A. if HP=42 or HP≥ 46)
- 20" 36 spokes 12" pneumatic*
if HP44; only curved rear frame)
- 24" MTB citybike tyre
- 24" MTB offroad tyre
- 24" Light (24 spokes)

- Pneum. Plugged (x 1"3/8) **Tyres**
 Pneum. high press. Plugged Smooth
- Solid x 1" 3/8
 - Solid KIK 24"
 - 24"/22" Pneum. Marathon Plus

Pushrims (14)

- Aluminium silver
 - Chromium-plated steel 24"
 - Aluminium dark
 - Rubber coated
 - OVAAL al.silver (*)
 - OVAALS al.dark (*)
 - OVAALG al.dark w/sil. insert (*)
 - Carbolife CURVE L 24"(*)
 - Carbolife GEKKO 24" (*)
 - Silicon covers 24"
- (*) not available with MTB wheels

Spoke guards (13)



- Transparent
 - Black
 - Fantasy
- code pattern:
 (see catalogue of spokeguards)

Front castors (14)

- 100 (4"x1.1/4) solid - plastic rim
- 100 (4"x1.1/4) solid - alum. rim
- 125 (5"x1.1/4) solid - plastic rim
- 125 (5"x1.1/4) solid - alum. rim
- 150 (6") pneumatic
- 150 (6") solid
- 175 (7") solid
- 175 (7") pneumatic

Straight Rear Frame

- With small wheels for narrow passages (default)
- Without any small wheels for narrow passages




Unbalancing system

- Right
 - Left
- (Only with straight rear frame)




Anti-tip device (15)

- Right
 - Left
- (Only with straight rear frame)




Single drive (16)

- (only with 24" or 22" rear wheels)
- Right
- Left




CURVED Rear Frame **Unbalancing system**



- Right
 - Left
- (Only with curved rear frame)

Revolving anti-tip device

- Right
 - Left
- (Only with curved rear frame)



(With curved rear frame it is not possible to assemble the small wheels for narrow passages)


Push handles (17)

- Standard Height adjustable
(max L = 48 cm)
(Compulsory with HP44cm or lower)




Centre of gravity (COG) (18)

- 2 cm (extended plate)
 - 3,5 cm
 - 5,5 cm
 - 7,5 cm
- (Measurements consider HP=44 cm)



Lateral Supports, swing away (19)

- Right
 - Left
- Specify the height "h" from the center of the swing away to the seat



Head-rest (20)

- Padded**
-
- Fixed
 - Height adjustable


- Orbital**
- Anatomic



- Moulded



- With lateral supports




Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = cm d = cm


Table (21)

- Transparent
- RxOLx Hemiplegic transp.
- Plastic grey
- Soft padded

- Specify the fastening:
- single double



Lined cushion in foam polyurethane (24)



- 3 cm high
- 5 cm high


Brakes (22)

- Classic - short handle
- Classic - long handle
- Drum
- Right single brake
- Left single brake
- Drum+classic

Options

- 25 "M" velcro strip (hook) on seat
- 26 "M" velcro for Vicair
- 27 "F" velcro strip (loop) on seat
- 28 Seat/back-rest connection flap
- 29 Chest harness
- 30 Chest belt
- 31 Abductor belt
- 32 Front protective covers
- 33 Stretched bar
- 34 Catheter holder pocket

Abductor (23)



- 1 3 5 7 cm
 - cm
- Choose or specify the height from seat upholstery

- 35 R L Crutch hold./single
- 36 R L Crutch hold./double
- 37 Pump for pneumatics
- 38 Suitcase holder
- 39 45° pelvic belt
- 40 45° padded pelvic belt
- 41 Black backpack
- 42 IV bar + poles
- 43 IV poles
- 44 Q.rel.axles Ø12mm w/42mm bushes (def.)
- 45 Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
- 46 Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes