

Order N° _____ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

Dealer _____

Key standard
 option

PAGE 1/2
 May 2024

L RANGE 1 W 28 30 32 34 36 38 40 42 44 46 48 50 52
 RANGE 2 W P > 36cm 32 34 36 38 40 42 44 46 48 50 52
 RANGE 3 W P > 36cm 32 34 36 38 40 42 44 46 48 50 52

P Seat depth Adjustable up to 6cm (max.50cm) 36 38 40 42 44 46
 48 50

HS Back-rest height 29 31 33 35 37 39 41 45
 47 49 51 55 60 (included stretched bar that can be enlarged)

HA Front seat/floor height with MAX size of front castors (fc) (NOTE: 175 front castors do not allow 82° detachable footrests)
 38(fc80max) 39(fc100max) 40(fc125max) 42(fc150max) 45(fc175max)
 44(fc150max) 46(fc175max)


HP Rear seat/floor height 34 36 38 40 42 44 45 46 48
 NOTE: 24" rear wheels (if HP min.38/max 48) 22" (if HP min.36/max 45) 20" (if HP min.34/max 43)

DP Foot-rest/seat distance 34 35 37 39 41 43 45
 (from 17 to 32 cm footrest directly assembled on the frame) cm

ALTHEA19PLUS

LOWERED SET-UP

Width & depth adjust.
 foldable alum.frame



Frame colour Matte black ①
 Metallic blue Pearl red Metallic grey
 Pearl white Glossy black
from: Neon colour
(see "Colours" catalogue)

Backrest upholstery ②
 Tension adj. nylon. back
 Tension adj. breath. back
 Without any backrest upholstery
 With divided seat (joint by velcro, def.)
Seat upholstery can grow +6cm max

Upholstery colour ③
 * (Coloured uphol. only with nylon)
 Black
 Red* Orange*
 Light blue* Green*
 Cornflower blue*

Not Detachable foot-rest ④

Detachable foot-rest 82° angle (DP 41 cm) (STANDARD) ④

Detachable foot-rest 68° angle (DP 41 cm) ④

Elevating foot-rest (detachable) ④

Plastic separate ⑤

Cover for footplates (recommended for seat width >52cm)
 Foot-holding belts (cover plates N.A.)
 Heel stops (cover plates N.A.)

Separate aluminium profile tilt adjustable footplates: ⑥

forward backward

Single aluminium profile tilt adjustable footplate: ⑦
(limited range of width adjustment)
 (not available with elevating foot-rest)

forward backward

Foot-holding belts
 Heel stops

Side-guards ⑧

Plastic Aluminium (same colour)

Detachable DESK armrest - tip-up ⑨

OPEN H22 with side-guard (only if HP ≥44 cm)
 OPEN H27 with side-guard (if HP ≤42cm)
 CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
 CLOSED H27 cm from seat (compulsory for HP=42cm&R24")

Elbow-rest
 Classic (23 cm long)
 Straight (26 cm long)
 Tubular (only if open armrest)
 Hemiplegic RX
 Hemiplegic LX
note: see also Hemi modular armrest ord.form

Detachable SPORT armrest - tip-up ⑩

OPEN H22 with side-guard (only if HP ≥44 cm)
 OPEN H27 with side-guard (if HP ≤42cm)
 CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
 CLOSED H27 cm from seat (compulsory for HP=42cm&R24")

Elbow-rest
 Classic (23 cm long)
 Straight (26 cm long)
 Tubular (only if open armrest)
 Hemiplegic RX
 Hemiplegic LX
note: see also Hemi mod.armrest ord.form

Detachable "U" armrest - tip-up ⑪

OPEN H22 with side-guard (only if HP ≥44 cm)
 OPEN H27 with side-guard (if HP ≤42cm)
 CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
 CLOSED H27 cm from seat (compulsory for HP=42cm&R24")

Elbow-rest
 Straight (26 cm long)
 Straight (36 cm long)
 Tubular (only if open armrest)
 Hemiplegic RX
 Hemiplegic LX
note: see also Hemi mod.armrest ord.form

Height adjustable tip-up arm-rest (only closed) ⑫

Desk model
 Sport model
 "U" model

Elbow-rest
 Classic (30 cm long)
 Classic (36 cm long)
 Straight (26 cm long)
 Straight (36 cm long)
 Hemiplegic RX
 Hemiplegic LX
note: see also Hemi mod.armrest ord.form

Dealer _____

Rear wheels (13)

24" 36 spokes 12" solid*
 22" 36 spokes 12" pneumatic*
 20" 36 spokes (N.A. if HP=42 if HP44: only curved rear frame)
 24" MTB citybike tyre
 24" MTB offroad tyre
 24" Light (24 spokes)

Pneu. Plugged (x 1"3/8) **Tyres**
 Pneu. high press. Plugged Smooth
 Solid grey x 1" 3/8
 Solid black KIK 24"
 24" Pneu. Marathon Plus

Pushrims (14)

Aluminium silver
 Chromium-plated steel 24"
 Aluminium dark
 Rubber coated
 OVAAL al.silver (*)
 OVAALS al.dark (*)
 OVAALG al.dark w/sil. insert (*)
 Carbolife CURVE L 24"(*)
 Carbolife GEKKO 24" (*)
 Silicon covers 24"
 (*) not available with MTB wheels

Spoke guards (15)

Transparent
 Black
 Fantasy
 code pattern:
 (see catalogue of spokeguards)

Front castors (16)

100 (4"x1.1/4) solid - plastic rim
 100 (4"x1.1/4) solid - alum. rim
 125 (5"x1.1/4) solid - plastic rim
 125 (5"x1.1/4) solid - alum. rim
 150 (6") pneumatic
 150 (6") solid
 175 (7") solid
 175 (7") pneumatic

Rear frame (17)

STRAIGHT rear frame with small wheels for narrow passages (default)
 STRAIGHT rear frame WITHOUT any small wheels f/narrow passages
 STRAIGHT rear frame with LEVER-activated narrow passage wheels
 CURVED rear frame

(Only if rear seat/floor height HP is between 44 and 48 cm)
 (It is not possible to assemble the small wheels for narrow passages)

Unbalancing system (18)

Right Left
 (Only with straight rear frame)
 Right Left
 (Only with curved rear frame)

Anti-tip device (Only with straight rear frame)
 Right
 Left

Revolving anti-tip device (Only with curved rear frame)
 Right
 Left

Single drive (19)

Right
 Left

(camber can be only 0° (rear wheels can be only 22" or 24") (the axle can grow +2cm or +4cm)

Extended wheel plate (20)

Push handles (21)

Standard Height adjustable Hinged backrest (no stretched bar) (max L=50)

Distance between the hub centre and back-rest (22)

-2 cm (ext. wheel plate)
 3,5 cm
 5,5 cm
 7,5 cm

(with HP=44 cm)

Adjustable swing away Lateral supports (23)

(tension adj. upholstery is compulsory)

Right
 Left

Specify the height "h" from the center of the swing away to the seat
 cm

Head-rest (24)

Padded (cannot grow in width)
 Fixed
 Height adjustable

Orbital
 Anatomic
 Moulded
 With lateral supports

Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = cm d = cm

Table (25)

Transparent
 RxOLx Hemiplegic transp.
 Plastic grey
 Soft padded

Specify the fastening:
 single double

(Only with arm-rest) if L>50cm

Lined cushion in foam polyurethane (28)

(cannot grow in width)

3 cm high
 5 cm high

Brakes (26)

Classic - short handle
 Classic - long handle
 Drum
 Right single brake
 Left single brake
 Drum+classic

Options

29 "M+F" velcro on seat
 30 "M" velcro strip (hook) on seat
 31 "F" velcro strip (loop) on seat
 32 Seat/back-rest connection flap
 33 Chest harness
 34 Chest belt
 35 Abductor belt
 36 Front protective covers
 37 Stretched bar (width adjust.)
 38 R OL Crutch hold./single
 39 R OL Crutch hold./double
 40 Pump for pneumatics

Abductor (27)

01 03 05 07 cm
 cm

Choose or specify the height from seat upholstery

41 45° pelvic belt
 42 45° padded pelvic belt
 43 Black backpack
 44 IV bar + poles
 45 IV poles

Quick rel. axles Ø 12 mm with 42mm cambered bushes (default)
 Quick rel. axles Ø 1/2" (12,7 mm) with 42mm cambered bushes
 Quick rel. axles Ø 1/2" (12,7 mm) with 50mm cambered bushes