

Order N° \_\_\_\_\_  Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer Ref. \_\_\_\_\_ User's weight Kg \_\_\_\_\_

Dealer \_\_\_\_\_

**Key**  standard  
 option


**PAGE 1/2**  
 May 2024

- L** Seat width 34 36 38 40 42 44 46 48 50  
52\* 56\* 60 (including eventual stretched bar, tension adj. upholstery & cover for footplate according to setup)
- P** Seat depth 36 38 40 42 44 46  
 Adjustable up to 6 cm (max.50cm)  
48 50
- HS** Back-rest height 29 31 33 35 37 39 41 43 45  
47 49 51 55 60 (included stretched bar)
- HA** Front seat/floor height 38(fc80max) 39(fc100max) 40(fc125max) 42(fc150max) 45(fc175max)  
 (with MAX size of front castors)  
 (NOTE: 175 front castors do not allow 82° detachable footrests)  
44(fc150max) 46(fc175max)
- HP** Rear seat/floor height 34 36 38 40 42 44 45 46 48  
 NOTE: 24" rear wheels (if HP min.38/max 48) 22" (if HP min.36/max 45) 20" (if HP min.34/max 43)
- DP** Foot-rest/seat distance 34 35 37 39 41 43 45  
 (from 17 to 32 cm footrest directly assembled on the frame) ..... cm

# ALTHEA.19

## LOWERED SET-UP

alum.foldable frame



- Frame colour**  Matte black ①  
 Metallic blue  Pearl red  Metallic grey  
 Pearl white  Glossy black  
 ..... Neon colour .....  
(see "Colours" catalogue)

- Upholstery** ②  
 Nylon (up to L50max)  
 Tension adj. nylon. backrest up.  
 Tension adj. breath. backrest up.  
 Without any backrest upholstery  
 With divided seat (joint by velcro, def.)  
 With one entire seat (not depth adjustable)

- Upholstery colour** ③  
 \* (Coloured uphol. only with nylon)  
 Black  
 Red\*  Orange \*  
 Light blue\*  Green \*  
 Cornflower blue \*

- Not Detachable foot-rest** ④
- Detachable foot-rest 82° angle** ④  
 (DP 41 cm)  
 (STANDARD)
- Detachable foot-rest 68° angle** ④  
 (DP 41 cm)
- Elevating foot-rest (detachable)** ④
- 

- Plastic separate** ⑤
- 
- Cover for footplates\*  
 Foot-holding belts (up to W50)  
 Heel stops (up to L50)  
 \* Eventual cover f/footplates included for W52-60

- Separate aluminium profile tilt adjustable foot-plates:** ⑥
- forward  backward
- 
- Foot-holding belts (up to W50)  
 \* Eventual cover f/footplates included for W52-60

- Single aluminium profile tilt adjustable foot-plate:** ⑦  
 (not available with elevating foot-rest)
- forward  backward
- 
- Foot-holding belts  
 Heel stops

- Automatic closure foot-plate in carbon** ⑧  
 (not available with detachable or elevating foot-rest; only L= from 36 to 46 cm)
- 
- Foot-holding belts

- Side guards** ⑨
- 
- Plastic
- Aluminium (same colour)

- Detachable DESK armrest - tip-up** ⑩
- OPEN H22 with side-guard (only if HP ≥44 cm)  
 OPEN H27 with side-guard (if HP≤42cm)
- CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
 CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
- Classic (23 cm long)  
 Straight (26 cm long)  
 Tubular (only if open armrest)  
 Hemiplegic RX  
 Hemiplegic LX  
 note: see also Hemi modular armrest ord.form

- Detachable SPORT armrest - tip-up** ⑪
- OPEN H22 with side-guard (only if HP ≥44 cm)  
 OPEN H27 with side-guard (if HP≤42cm)
- CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
 CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
- Classic (23 cm long)  
 Straight (26 cm long)  
 Tubular (only if open armrest)  
 Hemiplegic RX  
 Hemiplegic LX  
 note: see also Hemi mod.armrest ord.form

- Detachable "U" armrest - tip-up** ⑫
- OPEN H22 with side-guard (only if HP ≥44 cm)  
 OPEN H27 with side-guard (if HP≤42cm)
- CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
 CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
- Straight (26 cm long)  
 Straight (36 cm long)  
 Tubular (only if open armrest)  
 Hemiplegic RX  
 Hemiplegic LX  
 note: see also Hemi mod.armrest ord.form

- Height adjustable tip-up armrest (only closed)** ⑬
- Desk model  
 Sport model  
 "U" model
- Elbow-rest**
- Classic (30 cm long)  
 Classic (36 cm long)  
 Straight (26 cm long)  
 Straight (36 cm long)  
 Hemiplegic RX  
 Hemiplegic LX  
 note: see also Hemi modular armrest ord.form

Dealer \_\_\_\_\_

**Rear wheels** (14)

24" 36 spokes  12" solid\*  
 22" 36 spokes  12" pneumatic\*  
 20" 36 spokes (N.A. if HP=42 if HP44: only curved rear frame)  
 24" MTB citybike tyre  
 24" MTB offroad tyre  
 24" Light (24 spokes)

Pneum. Plugged (x 1"3/8) **Tyres**  
 Pneum. high press.  Plugged  Smooth  
 Solid grey x 1" 3/8  
 Solid black KIK 24"  
 24" Pneum. Marathon Plus

**Pushrims** (15)

Aluminium silver  
 Chromium-plated steel 24"  
 Aluminium dark  
 Rubber coated  
 Carbolife CURVE L 24"(\*)  
 OVAAL silver alum. (\*)  
 OVAALS dark alum. (\*)  
 OVAALG dark w/sil. insert (\*)  
 Carbolife GEKKO 24" (\*)  
 Silicon covers 24"  
 (\*) not available with MTB wheels

**Spoke guards** (16)

Transparent  
 Black  
 Fantasy  
 code pattern: .....  
 (see catalogue of spokeguards)

**Front castors** (17)

100 (4"x1.1/4) solid - plastic rim  
 100 (4"x1.1/4) solid - alum. rim  
 125 (5"x1.1/4) solid - plastic rim  
 125 (5"x1.1/4) solid - alum. rim  
 150 (6") pneumatic  
 150 (6") solid  
 175 (7") solid  
 175 (7") pneumatic

**Rear frame** (18)

STRAIGHT rear frame with small wheels for narrow passages (default)  
 STRAIGHT rear frame WITHOUT any small wheels f/narrow passages  
 STRAIGHT rear frame with LEVER-activated narrow passage wheels  
 CURVED rear frame

(Only if rear seat/floor height HP is between 44 and 48 cm)  
 (It is not possible to assemble the small wheels for narrow passages)

**Unbalancing system** (19)

Right  Right  
 Left  Left  
 (Only with straight rear frame) (Only with curved rear frame)

**Anti-tip device** (Only with straight rear frame)  
 Right  
 Left

**Revolving anti-tip device** (Only with curved rear frame)  
 Right  
 Left

**Single drive** (20)

Right  
 Left

**Extended wheel plate** (21)

**Push handles** (22)

Standard  Height adjustable(\*)  Hinged backrest(\*) (no stretched bar)

(\*) tension adj. upholstery is compulsory

**Distance between the hub centre and back-rest** (23)

-2 cm (ext. wheel plate)  
 3,5 cm  
 5,5 cm  
 7,5 cm

(with HP=44 cm)

**Adjustable swing away Lateral supports** (23)

(tension adj. upholstery is compulsory)

Right  
 Left

Specify the height "h" from the center of the swing away to the seat  
 ..... cm

**Head-rest** (24)

**Padded**

Fixed  
 Height adjustable

**Orbital**

Anatomic  
 Moulded  
 With lateral supports

Specify "h" & "d" where requested  
 Standard: h = 20 cm d = 15 cm  
 h = ..... cm d = ..... cm

**Table** (25)

Transparent  
 RxOLx Hemiplegic transp.  
 Plastic grey  
 Soft padded

Specify the fastening:  
 single  double

(Only with arm-rest)

**Lined cushion in foam polyurethane** (28)

3 cm high  
 5 cm high

**Brakes** (26)

Classic - short handle  
 Classic - long handle  
 Drum  
 Right single brake  
 Left single brake  
 Drum+classic

**Options**

29  "M" velcro strip (hook) on seat  
 30  "M" velcro strip for Vicair  
 31  "F" velcro strip (loop) on seat  
 32  Seat/back-rest connection flap  
 33  Chest harness  
 34  Chest belt  
 35  Abductor belt  
 36  Front protective covers  
 37  Stretched bar  
 38  Catheter holder pocket  
 39  R  L Crutch hold./single  
 40  R  L Crutch hold./double

**Abductor** (27)

1  3  5  7 cm  
 ..... cm  
 Choose or specify the height from seat cover

41  Pump for pneumatics  
 42  Suitcase holder  
 43  45° pelvic belt  
 44  45° padded pelvic belt  
 45  Black backpack  
 46  IV bar + poles  
 47  IV poles  
 48  Padded rigid cushion (L52max)  
 49  Padded rigid cushion (L56-60)  
 50  Q.rel.axles Ø12mm w/42mm bushes (def.)  
 51  Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes  
 52  Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes  
 53  **Stickers for car fastenings**  
 (escl.: load>120Kg, R12", UltraL.RW, hinged bckr)