


Dealer

Key standard option

ANTARES

- L Seat width 30 32 34 36 38 40 42 44
46 48 50 (double gas spring for tilt-in & reclining blackrest)
- P Seat depth 34 36 40 42 44 46
- HS Back-rest height 37 41 45 50 55
- DP Foot-rest/seat distance 31 33 35 37 39 41 43 45 47 49 cm



- Frame colour** Metallic blue Metallic grey (see "Colours" catalogue) neon colours
- Upholstery colour** (only nylon) Black Red Light blue Cornflower blue
- Pearl red Matte black Orange Green

- Seat upholstery** Without seat Nylon upholstery Rigid base without upholstery Anatomic seat in microfiber, up to L 46 cm Anatomic seat in microfiber, L > 46 cm

- Backrest upholstery** Without backrest Tension adjustable nylon upholstery Tension adjustable breathable upholstery Anatomic backrest in microfiber, up to L 46 cm Anatomic backrest in microfiber, L > 46 cm

- Foot-rest** Fixed Detachable Detachable with spring Detachable elevating
-

- Plastic separate** (only L from 38 to 46 cm) (not available with foot-rest with spring) Foot-holding belts Heel stops
-

- Separate aluminium profile footplates** (not available with foot-rest with spring) Foot-holding belts
-

- Height adjustable armrest** Specify the height "h" of the arm-rest from the seat
-

- Single aluminium profile footplate** (not available with elevating and spring foot-rest) Foot-holding belts with locking mechanism
-

- Brakes** Classic - short handle Classic - long handle Drum Classic + Drum

- Rear wheels** 12" (300mm) pneumatic 12" (300mm) solid 20" 22" 24"
- Tyres** Pneumatic plugged Solid 20"/22" (x 1.3/8") Solid KIK 24" black

- Handrims** (only upon request) Not available with 12" wheels Aluminium Chromium-plated steel Rubber coated Silicon covers

- Spoke guards** Not available with 12" wheels Transparent Black Fantasy code pattern: (see catalogue of spokeguards)

- Front castors** 125 (5") solid 150 (6") solid 175 (7") pneumatic 175 (7") solid

