

☐ Order N ° _____ ☐ Quotation
 Date _____ Signature _____
 Customer ref. _____

PAGE 1/2
 Jan 2024

Dealer _____


Key

☐ standard

☐ option

- L** Seat width ☐ 24 ☐ 26 ☐ 28 ☐ 30 ☐ 32 ☐ 34 ☐ 36
- P** Seat depth ☐ 26 ☐ 28 ☐ 30 ☐ 32 ☐ 34 ☐ 36
- HS** Back-rest height ☐ 21 ☐ 23 ☐ 25 ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39
- HA** Front seat/floor height with MAX size of front castors (fc) ☐ 38(fc80max) ☐ 39(fc100max) ☐ 40(fc125max) ☐ 42(fc150max) ☐ 45(fc175max) ☐ 48
☐ 44(fc150max) ☐ 46(fc175max)
- HP** Rear seat/floor height ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
 NOTE: 24" rear wheels if min.HP 38-max HP 48 22" if min.HP 36- max HP 45 20" if min.HP 34- max HP 43
- DP** Foot-rest/seat distance ☐ 25 ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 cm

ARIEL



Frame colour

①

- ☐ Met. blue ☐ Pearl red ☐ Pastel yellow
☐ Deep sky light blue ☐ Met. pink ☐ Matte black
☐ Neon orange ☐ Neon green
☐ (see "Colours" catalogue)

Upholstery

②

- ☐ Tension adj. nylon back-rest
☐ Tension adj. breathable backr.
☐ Without back-rest
☐ Without seat

Uphostery colour

③

- (nylon only)
☐ Black ☐ Light blue
☐ Red ☐ Green
☐ Cornflower blue ☐ Orange

Footrests


④

- ☐ Not detachable footrests, **STRAIGHT** front frame (Compulsory pull to brake if P<30cm)
☐ Detachable footrests, **STRAIGHT** front frame (Compulsory pull to brake)
☐ **ABDUCTED** front frame with NOT detachable foot-rests (Need technical approval)
☐ **ABDUCTED** front frame with detachable foot-rests (Need technical approval)
☐ **Elevating footrests** (Need technical approval)
 (Std brake starting from P32cm)

☐ **Separate aluminium profile foot-plates**

⑤

(from DP 17 to DP 31 cm the footplate is assembled on the frame)
 (for L24-26-28cm recommended single footplate)




☐ Foot-holding belts

☐ **Single aluminium profile tilt adjust. footplate**

⑥

(from DP 17 to DP 31 cm the footplate is assembled on the frame)



☐ Foot-holding belts

Rear wheel camber


⑦

- ☐ 0° ☐ 1° ☐ 3° ☐ 5°

- ☐ Quick rel. axles Ø 12 mm (default)
☐ Quick rel. axles Ø 1/2" (12,7 mm)
☐ Cambered bushes length 42 mm (def.)
☐ Cambered bushes length 50 mm

Side guards

☐ Plastic



☐ Aluminium (same colour)

☐ **Height adjustable "T" armrest (with plastic side-guards)**

⑧




Elbow-rest

- ☐ Classic (23 cm long)
☐ Classic (30 cm long)
☐ Straight (26 cm long)

☐ **Height adjustable "T" armrest (with aluminium side-guards)**

⑨




Elbow-rest

- ☐ Classic (23 cm long)
☐ Classic (30 cm long)
☐ Straight (26 cm long)

☐ **Tip-up height adjustable from 22 to 28cm "L" armrest**

⑩



Elbow-rest

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Rear wheels

⑫

- ☐ 24" (36 spokes) (*)
☐ 22" (36 spokes)
☐ 20" (36 spokes)
☐ 24" Light (24 spokes) (*)
☐ 22" Light (24 spokes)
☐ 20" Light (20 spokes)
☐ 12" (300) pneumatic (**)
☐ 12" (300) solid (**)
 (**) need technical approval
Tyres
☐ Pneum. HP smooth 20"/22"/24"
☐ Pneum. HP plugged 22"/24"
☐ Solid 24x1.3/8"
☐ Solid KIK 24" black
☐ Marathon Plus 22" or 24"

(*) Only with P=36 or P34cm

Pushrims

⑬

- ☐ Aluminium silver
☐ Chromium-plated steel 24"
☐ Aluminium dark
☐ Rubber coated
☐ OVAAL alum. silver 24" (*)
☐ OVAALS alum. dark 24" (*)
☐ OVAALG dark w/silic. insert 24" (*)
☐ Carbolife CURVE L 24" (*)
☐ Carbolife GEKKO 24" (*)
☐ Silicon covers

(*) not available with MTB wheels

Spoke guards

⑭

- ☐ Transparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)
☐ Fantasy with name
 code pattern:
 (see catalogue of spokeguards)
 name to be printed:
 "....."

Front castors

⑮

- ☐ 80 (3") solid with lights
☐ 80 (3") solid
☐ 100 (4") solid with lights
☐ 100 (4") solid
☐ 100 (4") solid with alum. rim
☐ 125 (5") solid
☐ 125 (5") solid with lights
☐ 125 (5") solid with alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 Jan 2024


ARIEL

Dealer


Push handles

17


☐ Standard



☐ Height adjustable




☐ Abducted height adjustable




☐ Curved rear frame


(it is not possible to assemble the small wheels for narrow passages)



☐ With straight rear frame



☐ With small wheels for narrow passages



(need technical approval)

Single drive

(camber can be only 0°)
 (rear wheels can be only 22" or 24")


21



☐ Right
☐ Left

☐ Extended wheel plate

22



Unbalancing system

20



☐ Right

☐ Left




☐ Right

☐ Left

(Only with curved rear frame)

Anti-tip device



☐ Right

☐ Left

☐ Right & left

Revolving anti-tip device

(Only with curved rear frame)



☐ Right


☐ Left

☐ Right & left

Head-rest

23

Padded




☐ Height adjustable

Orbital


☐ Anatomic



☐ Moulded




☐ With lateral supports



Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = cm d = cm

☐ Abductor

24



☐ 1 ☐ 3 ☐ 5 ☐ 7 cm

☐ cm

Choose or specify the height from the seat cover

Anodized set

(Rear wheel rims & hub if Light and forks)

27

☐ Black

☐ Silver

☐ Light blue

☐ Red

☐ Orange

☐ Cyclamen light purple

(Light wheels only)

Table

26

☐ Transparent

☐ Soft padded


Specify the fastening:

☐ single ☐ double

(Only with arm-rest)

Adjustable lateral supports

25



☐ Right

☐ Left

Specify the height "h" from seat cover to the center of the swing away

Lined cushion

in foam polyurethane

28



☐ 3 cm high

☐ 5 cm high

Brakes

29

☐ Classic - short handle (from P32cm)

☐ Classic - pull to brake

☐ Classic long handle

☐ Brakes inside the sideguard*

☐ Right single brake

☐ Left single brake

☐ Drum

☐ Drum + classic

* compulsory for some configurations

Options

- 30 ☐ "M" velcro (hook) on seat
 31 ☐ "M" velcro for VICAIR cushion
 32 ☐ "F" (loop) on seat
 33 ☐ Seat/backrest connection flap
 34 ☐ Padded rigid cushion
 35 ☐ Chest harness
 36 ☐ Chest belt
 37 ☐ Abductor belt
 38 ☐ Double cross bar *
 39 ☐ Stretched bar
 40 ☐ Pump for pneumatics
 41 ☐ 45° pelvic belt
 42 ☐ 45° padded pelvic belt
 43 ☐ Black backpack
 44 ☐ Front protective covers

* need technical approval

Notes