

Dealer _____

Key ☐ standard ☐ option☐ **ARYA** (seat depth adjustable only)☐ **ARYA PLUS** (seat width & depth adjustable)

L Seat width ☐20 ☐22 ☐24 ☐26
☐24 ☐26 ☐28 ☐30 ☐32
☐28 ☐30 ☐32 ☐34 ☐36

P Seat depth ☐20 ☐22 ☐24 ☐26 ☐28
 (adjustable +2, +4, +6 cm ☐30 ☐32 ☐34 ☐36 ☐38
 from any chosen measurement)


HS Back-rest height ☐25 ☐27 ☐29 ☐31 ☐33 ☐35 ☐37 ☐39 ☐41

HA Front seat/floor height ☐40 ☐42 ☐44 ☐46 ☐48

HP Rear-seat/floor distance ☐36 ☐38 ☐40 ☐42 ☐44 ☐46

DP Foot-rest/seat distance cm (min. 15 cm, max 35 cm for HA=40 cm)
 (min. 15 cm, max 41 cm for HA=46 cm)

cm

ARYA**Frame colour**

- ☐ Pearl red ☐ Pearl white ☐ Metallic pink
☐ Metallic light blue ☐ Orange FLUO ☐ Green FLUO
☐ (see "Colours" catalogue)
☐ Two-colour version a) PEARL WHITE sitting frame
 b) colour of the adjustment bar: (see "Colours" catalogue)

①

Backrest upholstery

- ☐ Tension adjustable nylon
☐ Tension adj. breath. back
☐ Without back-rest

②

Upholstery colour

(Nylon only)

- ☐ Black ☐ Light blue
☐ Red ☐ Green
☐ Conflower blue ☐ Orange

③

Front frame

④


☐ Standard☐ Standard
(not suggested
for L<32)☐ Abducted**Foot-rest**

⑤

- ☐ Not detachable
foot-rest
☐ Detachable
foot-rest

☐ **Single aluminium
footplate**


⑥

☐ Foot-holding belts☐ **Separate aluminium
footplate**

Only for L>32 cm (standard front frame)
 Only for L>28 cm (abducted front frame)

☐ Foot-holding belts

⑦





☐ Open

**Liftable single
height & depth adjustable
foot-plate**

Single choice of footplate for front open frame

☐ Standard
push
handles☐ Straight height
adjustable
push
handles)☐ Abducted height
adjustable
push
handles☐ Height
adjustable
push
bar

⑧

**Side-guards**☐ Plastic☐ Carbon fibre☐ Aluminium
(same colour)

⑨


☐ **Height adjustable
"T" integrated armrest**
(to add to side-guards)**Elbow-rest**

- ☐ Classic (23 cm long)
☐ Classic (30 cm long)
☐ Straight (23 cm long)

⑩

**Distance between the hub
centre and backrest**


⑪



- ☐ -3 cm (extended wheel plate)
☐ 0 cm
☐ 3 cm (standard)
☐ 5 cm

Backrest tilt


⑫



- ☐ 89°
☐ 91° (standard)
☐ 93°

Revolving anti-tip device

⑬



- ☐ Right
☐ Left

Note

Dealer _____

Rear wheels

⑭

- ☐
- 24" Light(24 spokes) (*)
-
- ☐
- 22" Light(24 spokes)
-
- ☐
- 20" Light(20 spokes)

Tyres

- ☐
- Pneum. high pressure smooth
-
- ☐
- Pneum. high pressure plugged
-
- ☐
- Pneum. plugged (24"x1.3/8")
-
- ☐
- Solid 24x1.3/8"
-
- ☐
- Solid KIK 24" black
-
- ☐
- Marathon Plus 22" or 24"

 (*) Only with P=36 cm and
 HP from 40 to 46 cm

Rear wheel camber

⑮

- ☐
- 1°
- ☐
- 3°
- ☐
- 5°

Pushrims

⑮

- ☐
- Aluminium silver
-
- ☐
- Chromium-plated steel 24"
-
- ☐
- Dark aluminium
-
- ☐
- Rubber coated
-
- ☐
- Carbolife Curve L 24"
-
- ☐
- Carbolife GEKKO 24"
-
- ☐
- Silicon covers

(*) Not available with 20" wheels

Spoke guards

⑮

- ☐
- Trasparent
-
- ☐
- Black
-
- ☐
- Fantasy

 code pattern:
 (see catalogue of spokeguards)

- ☐
- Fantasy with name

 code pattern:
 (see catalogue of spokeguards)
 name to be printed:

"....."

Front castors


⑮

- ☐
- 80 (3") solid with lights
-
- ☐
- 80 (3") solid
-
- ☐
- 100 (4") solid with lights
-
- ☐
- 100 (4") solid
-
- ☐
- 100 (4") solid with alum. rim
-
- ☐
- 125 (5") solid
-
- ☐
- 125 (5") solid with alum.rim
-
- ☐
- 150 (6") pneumatic
-
- ☐
- 150 (6") solid

Head-rest

⑮

Padded



- ☐
- Height adjustable

Orbital

- ☐
- Anatomic


- ☐
- Moulded

- ☐
- With lateral supports


 Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = cm d = cm

Abductor

⑮




- ☐
- 1
- ☐
- 3
- ☐
- 5
- ☐
- 7 cm
- ☐
- cm

Pls choose or specify height from seat

Lateral supports

⑮

Swing Away



- ☐
- Right
-
- ☐
- Left

Height adjustable

(Pls specify h)

Table

⑮

- ☐
- Trasparent
-
- ☐
- Soft padded

Pls specify the fastening

- ☐
- single
- ☐
- double



(only with arm-rest)

Anodized set

⑮

(Rear wheel rims&hubs and forks)

- ☐
- Black
-
- ☐
- Silver
-
- ☐
- Light blue
-
- ☐
- Red
-
- ☐
- Orange
-
- ☐
- Cyclamen light purple

Lined cushion

⑮

in foam polyurethane



- ☐
- 3 cm high
-
- ☐
- 5 cm high

Brakes

⑮

- ☐
- Classic (inside the sideguard)
-
- ☐
- Pulling brake short handle inside the sideguard
-
- ☐
- Drum
-
- ☐
- Drum + classic (inside the sideguard)

Options

- 26
- ☐
- "M" velcro strip (hook) on seat
-
- 27
- ☐
- "F" velcro strip (loop) on seat
-
- 28
- ☐
- Padded rigid seat
-
- 29
- ☐
- Chest harness
-
- 30
- ☐
- Chest belt
-
- 31
- ☐
- Abductor belt
-
- 32
- ☐
- Pump for pneumatics
-
- 33
- ☐
- 45° pelvic belt
-
- 34
- ☐
- 45° padded pelvic belt
-
- 35
- ☐
- Black holdall
-
- 36
- ☐
- Front protective covers

Note