

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 1/2
 Jan 2024

Dealer _____

Key ☐ standard ☐ option

CHILDREN 3000

L Seat width

☐ 24* ☐ 26* ☐ 28* ☐ 30 ☐ 32 ☐ 34 ☐ 36

* Not available with detachable foot-rests

CHILDREN 3000 PLUS

L Seat width (adjustable within the selected range)

☐ RANGE 1 ☐ 24* ☐ 26* ☐ 28* ☐ 30 ☐ 32 ☐ 34

☐ RANGE 2 ☐ 26* ☐ 28* ☐ 30 ☐ 32 ☐ 34 ☐ 36 ☐ 38

* Not available with detachable foot-rests

P Seat depth

☐ 28 ☐ 30 (adjustable from 28 to 34 cm)

☐ 32 ☐ 34 ☐ 36 ☐ 38 (adjustable from 32 to 38 cm)

HS Back-rest height

☐ 25 ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41

HA Front seat/floor height

☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48

HP Rear seat/floor height

☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46

DP Foot-rest/seat distance

☐ 17* ☐ 19* ☐ 21* ☐ 23 ☐ 25 ☐ 27 ☐ 29 ☐ 31 ☐ 33


☐ 35 ☐ 37 ☐ 39

* Not available with detachable foot-rests (unless the front frame is abducted)

cm


CHILDREN 3000

Fold.alum.frame for Kid



☐ with CHILDREN 3000 rear frame

☐ with CURVED rear frame



Frame colour

- ☐ Deep sky light blue ☐ Orange ☐ Pearl red ☐ Metallic pink
☐ Neon orange ☐ Neon green
☐ (see "Colour" catalogue)

①

Backrest upholstery

- ☐ Tension adjustable nylon back-rest
☐ Tension adj. breathable b.
☐ Without back-rest
☐ Without seat

②

Upholstery colour (only nylon)


- ☐ Black ☐ Light blue
☐ Red ☐ Green
☐ Cornflower blue ☐ Orange

③


Front frame

④

☐ Standard
 (Not suggested for L < 32cm)



☐ Abducted



Foot-rest

⑤


☐ Detachable

☐ Not detachable

☐ Elevating (detachable)
 Subject to technical approval

Separate aluminium profile footplates


⑥



☐ Foot-holding belts

Single aluminium profile footplate


⑦



☐ Foot-holding belts

Side-guards

⑧



☐ Plastic (standard)


(Not available with reclining backrest)



☐ Aluminium (curved, same colour)

Height adjustable "T" armrest

⑨



Elbow-rest




- ☐ Classic (23 cm long)
☐ Classic (30 cm long)
☐ Straight (26 cm long)

Standard push handles

Straight height adjustable push handles

Abducted height adjustable push handles

⑩



Rear wheels

⑪

- ☐ 24" Light (24 spokes) (*)
☐ 22" Light (24 spokes)
☐ 20" Light (20 spokes)

Tyres

- ☐ Pneum. high pressure smooth
☐ Pneum. high pressure plugged
☐ Solid Grey 20"/22"x1/24x1.3/8
☐ Solid black KIK 24" x1
☐ Marathon Plus 22" or 24"
 (*) Only with P = 34 or 36 cm and HP from 40 to 46 cm

Handrims

⑫

- ☐ Aluminium silver
☐ Chromium-plated steel
☐ Dark aluminium
☐ Rubber coated
☐ Carbolife CURVE L 24"
☐ Carbolife GEKKO 24"
☐ Silicon covers

Spoke guards

⑬

- ☐ Transparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)

Front castors

⑭


- ☐ 80 (3") solid with lights
☐ 80 (3") solid
☐ 100 (4") solid with lights
☐ 100 (4") solid
☐ 100 (4") solid with alum. rim
☐ 125 (5") solid
☐ 125 (5") solid with lights
☐ 125 (5") solid with alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid

Rear wheel camber

- ☐ 1° ☐ 3° ☐ 5°

- ☐ Quick release axles Ø12mm with 42mm bushes (def.)
☐ Quick release axles Ø1/2" (12,7mm) with 42mm bushes
☐ Quick release axles Ø1/2" (12,7mm) with 50mm bushes


⑮

Liftable anti-tip device Children 3000
 for Children 3000 rear frame



- ☐ Right
☐ Left
☐ Right and left

Revolving anti-tip device
 for CURVED rear frame


- ☐ Right
☐ Left
☐ Right and left € 132.00

☐ **Extended wheel plate**



Rear frame of Children 3000 (requires liftable anti-tip for Children 3000 and is compulsory with reclining back-rest)



Curved rear frame (requires revolving anti-tip device)

☐ **Single drive**


(Only with 0° rear wheel camber)
 (Only with 22" or 24" wheels)



- ☐ Right
☐ Left

Head-rest

Padded




- ☐ Height adjustable

Orbital


- ☐ Anatomic



- ☐ Moulded




- ☐ With lateral supports




Specify "h" & "d" where requested
 Standard: h = 15 cm d = 10 cm

h = cm d = cm

Abductor


- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm
☐ cm

Choose or specify the height from the seat cover

Lateral supports


- ☐ Right
☐ Left

Specify the height "h"

Table

- ☐ Transparent
☐ Black soft padded

specify the fastening:

- ☐ single ☐ double



(Only with arm-rest)

Anodized set

(Rear wheel rim&hub and front fork)

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

Brakes

- ☐ Classic - short handle
☐ Classic - pulling brake
☐ Classic - long handle
☐ Right single brake
☐ Left single brake
☐ Drum
☐ Drum + classic

Lined cushion

in foam polyurethane



- ☐ 3 cm high
☐ 5 cm high


Options

- 28 ☐ "M" & "F" velcro on seat
 29 ☐ "M" velcro strip (hook) on seat
 30 ☐ "F" velcro strip (loop) on seat
 31 ☐ Label of connection seat/back-rest
 32 ☐ Padded rigid cushion
 33 ☐ Chest harness
 34 ☐ Chest belt
 35 ☐ Abductor belt
 36 ☐ Pump for pneumatics
 37 ☐ 45° pelvic belt
 38 ☐ 45° padded pelvic belt
 39 ☐ Black backpack
 40 ☐ **Stickers for car fastenings** **crash tested**

☐ **Reclining backrest**

from 90° to 150°
 (compulsory:

- rear frame for Children 3000
 - straight or abducted height adjustable push handles
 - stretched bar
 - extended wheel plate


Notes