


Dealer _____

- L** Seat width (adjustable within the selected range) RANGE 1 W 24* 26* 28* 30* 32
 RANGE 2 W 28* 30* 32 34 36
 * Recommended abducted front frame
- P** Seat depth 26 28 30 32 34 36
- HS** Back-rest height 21 23 25 27 29 31 33 35 37 39
- HA** Front seat/floor height with MAX size of front castors (fc) 38(fc80max) 39(fc100max) 40(fc125max) 42(fc150max) 45 48
 44(fc150max) 46
- HP** Rear seat/floor height NOTE: 24" rear wheels if min.HP 38-max HP 48 34 36 38 40 42 44 46 48
 22" if min.HP 36- max HP 45 20" if min.HP 34- max HP 43
- DP** Foot-rest/seat distance 25 27 29 31 33 35 37 39 41 cm

Key
 standard
 option


ARIEL PLUS

Width & depth adjust.
 KID fold.alum.frame




- Frame colour** Met. blue Pearl red Pastel yellow
 Deep sky light blue Met. pink Matte black
 Neon orange Neon green
 (see "Colours" catalogue)
- Backrest upholstery** Tension adj. nylon back-rest
 Tension adj. breathable backr.
 Without back-rest
- Uphostery colour** (Nylon only) Black Light blue
 Red Green
 Cornflower blue Orange

- Straight front frame** (Recommended for seat width L > 32 cm) Not detachable footrests
 Detachable footrests (Compulsory pull to brake if depth < 30cm)
 Elevating (Need technical approval) (Push brake starting from P32cm)




- Separate aluminium profile footplates** (from DP 17 to DP 31 cm the footplate is assembled on the frame) (for L24-26-28cm recommended single footplate)



Foot-holding belts

- Single aluminium profile tilt adjust. footplate** (from DP 17 to DP 31 cm the footplate is assembled on the frame)



Foot-holding belts

- Rear wheel camber** 0° 1° 3° 5°
- Q.rel.axles Ø12mm w/42mm bushes (def.)
 Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
 Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes

- Abducted front frame** (Recommended for seat width L < 32 cm) Abducted front frame with NOT detach.foot-rests
 Abducted front frame with detach.foot-rests
 Elevating (To be added to abducted frame) (Push brake starting from P32cm) (Need technical approval)



- Side-guards and ev. integrated armrests** Plastic side-guard Aluminium side-guard (frame colour)
 Height adjustable "T" armrest (with plastic side-guards) Height adjustable "T" armrest (with aluminium side-guards)



Elbow-rest

- Classic (23 cm long)
 Classic (30 cm long)
 Straight (26 cm long)

- Tip-up height adjustable from 22 to 28cm "L" armrest** (sideguard is needed)



Elbow-rest

- Tubular
 Straight (26 cm long)
 Classic (23 cm long)
 Classic (30 cm long)

- Rear wheels** 24" (36 spokes) (*) (* Only with P=36 or P34cm)
 22" (36 spokes)
 20" (36 spokes)
 24" Light (24 spokes) (**)
 22" Light (24 spokes)
 20" Light (20 spokes)
 12" (300) pneumatic (**)
 12" (300) solid (**)
 (**) need technical approval

Tyres

- Pneum. HP smooth 20"/22"/24"
 Pneum. HP plugged 22"/24"
 Solid Grey 20"/22"x1/24"x1.3/8"
 Solid Black KIK 24"
 Marathon Plus 22" or 24"

- Pushrims** Aluminium silver
 Chromium-plated steel 24"
 Aluminium dark
 Rubber coated
 OVAAL al.silver 24"
 OVAALS al.dark 24"
 OVAALG al.dark w/sil. insert 24"
 Carbolife CURVE L 24"
 Carbolife GEKKO 24"
 Silicon covers 24"


- Spoke guards** Trasparent
 Black
 Fantasy
 code pattern:
 (see catalogue of spokeguards)
 Fantasy with name
 code pattern:
 (see catalogue of spokeguards)
 name to be printed:
 "....."

- Front castors** 80 (3") solid with lights
 80 (3") solid
 100 (4") solid with lights
 100 (4") solid
 100 (4") solid with alum. rim
 125 (5") solid
 125 (5") solid with lights
 125 (5") solid with alum. rim
 150 (6") pneumatic
 150 (6") solid


Dealer

Push handles (15)


Standard



Height adjustable




Abducted height adjustable




Curved rear frame


(it is not possible to assemble the small wheels for narrow passages)



With straight rear frame





With small wheels for narrow passages




(need technical approval)

Unbalancing system (17)



Revolving anti-tip device

(Only with curved rear frame)




Anti-tip device

(Only with straight rear frame)




Single drive (18)

(camber can be only 0°)
 (rear wheels can be only 22" or 24")
 (the axle can grow +2cm or +4cm)




Right
 Left

Extended wheel plate (19)



Head-rest (20)

Padded
 (cannot grow in width)



Height adjustable

Abductor (21)

onto the frame for P > 28cm




onto the seat for P < 28cm

1 3 5 7 cm
 cm

Choose or specify the height from the seat cover

Adjustable lateral supports (22)



Right
 Left

Specify the height "h" from the center of the swing away to the seat

Orbital


Anatomic



Moulded



With lateral supports



Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = cm d = cm

Anodized set (23)

(Rear wheel rims & hub if Light and forks)

- Black
- Silver
- Light blue
- Red
- Orange
- Cyclamen light purple
 (Light wheels only)


Table (25)

- Transparent
- Soft padded

(Only with armrests and always with double fastening)

Lined cushion (24)

in foam polyurethane
 (cannot grow in width)



- 3 cm high
- 5 cm high

Brakes (26)

- Classic push (short handle) from P32cm
- Classic pull (short handle)
- Classic push long handle
- Classic pull long handle
- Brakes inside the sideguard*
- Right single brake
- Left single brake
- Drum
- Drum + classic

* compulsory for some configurations

Options

- 27 "M+F" velcro on seat
- 28 "M" velcro (hook) on seat
- 29 "M" velcro for VICAIR cushion
- 30 "F" (loop) on seat
- 31 Seat/backrest connection flap
- 32 Chest harness
- 33 Chest belt
- 34 Abductor belt
- 35 Stretched bar (width adj.)
- 36 Pump for pneumatics
- 37 45° pelvic belt
- 38 45° padded pelvic belt
- 39 Black backpack
- 40 Front protective covers

Notes