


Dealer _____

ARIEL



Key	<input type="checkbox"/> standard	<input type="radio"/> option
L Seat width	<input type="checkbox"/> 24 <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36	
P Seat depth	<input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36	
HS Back-rest height	<input type="checkbox"/> 21 <input type="checkbox"/> 23 <input type="checkbox"/> 25 <input type="checkbox"/> 27 <input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39	
HA Front seat/floor height with MAX size of front castors (fc)	<input type="checkbox"/> 38(fc80max) <input type="checkbox"/> 39(fc100max) <input type="checkbox"/> 40(fc125max) <input type="checkbox"/> 42(fc150max) <input type="checkbox"/> 45(fc175max) <input type="checkbox"/> 48 <input type="checkbox"/> 44(fc150max) <input type="checkbox"/> 46(fc175max)	
HP Rear seat/floor height NOTE: 24" rear wheels if min.HP 38-max HP 48 22" if min.HP 36-max HP 45 20" if min.HP 34-max HP 43	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48	
DP Foot-rest/seat distance	<input type="checkbox"/> 25 <input type="checkbox"/> 27 <input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 <input type="checkbox"/> 43	cm

Frame colour ①

Met. blue Pearl red Pastel yellow
 Deep sky light blue Met. pink Matte black
 Neon orange Neon green
 (see "Colours" catalogue)

Upholstery ②

Tension adj. nylon back-rest
 Tension adj. breathable backr.
 Without back-rest
 Without seat

Upholstery colour ③
(nylon only)

Black Light blue
 Red Green
 Cornflower blue Orange

Footrests ④

Not detachable footrests, **STRAIGHT** front frame
(Compulsory pull to brake if P<30cm)
 Detachable footrests, **STRAIGHT** front frame
(Compulsory pull to brake)
 ABDUCTED front frame with NOT detachable foot-rests
(Need technical approval)
 ABDUCTED front frame with detachable foot-rests
(Need technical approval)
 Elevating footrests
(Need technical approval)
(Std brake starting from P32cm)

Separate aluminium profile foot-plates ⑤
(from DP 17 to DP 31 cm the footplate is assembled on the frame)
(for L24-26-28cm recommended single footplate)
 Foot-holding belts

Single aluminium profile tilt adjust. footplate ⑥
(from DP 17 to DP 31 cm the footplate is assembled on the frame)
 Foot-holding belts

Rear wheel camber ⑦

0° 1° 3° 5°

Quick rel. axles Ø 12 mm (default)
 Quick rel. axles Ø 1/2" (12,7 mm)
 Cambered bushes length 42 mm (def.)
 Cambered bushes length 50 mm

Side guards

Plastic
 Aluminium (same colour)

Height adjustable "T" armrest (with plastic side-guards)

Elbow-rest

Classic (23 cm long)
 Classic (30 cm long)
 Straight (26 cm long)

Height adjustable "T" armrest (with aluminium side-guards) ⑨

Elbow-rest

Classic (23 cm long)
 Classic (30 cm long)
 Straight (26 cm long)

Tip-up height adjustable from 22 to 28cm "L" armrest ⑩

Elbow-rest

Tubular
 Straight (26 cm long)
 Classic (23 cm long)
 Classic (30 cm long)

Rear wheels ⑫

24" (36 spokes) (*) (* Only with P=36 or P34cm)
 22" (36 spokes)
 20" (36 spokes)
 24" Light (24 spokes) (*)
 22" Light (24 spokes)
 20" Light (20 spokes)
 12" (300) pneumatic (**)
 12" (300) solid (**)
 (**) need technical approval

Tyres

Pneum. HP smooth 20"/22"/24"
 Pneum. HP plugged 22"/24"
 Solid 24x1.3/8"
 Solid KIK 24" black
 Marathon Plus 22" or 24"

Pushrims ⑬

Aluminium silver
 Chromium-plated steel 24"
 Aluminium dark
 Rubber coated
 OVAAL alum. silver 24" (*)
 OVAALS alum. dark 24" (*)
 OVAALG dark w/silic. insert 24" (*)
 Carbolife CURVE L 24" (*)
 Carbolife GEKKO 24" (*)
 Silicon covers

(*) not available with MTB wheels

Spoke guards ⑭

Transparent
 Black
 Fantasy
 Fantasy with name
 code pattern:
 (see catalogue of spokeguards)
 name to be printed:
 "....."


Front castors ⑮

80 (3") solid with lights
 80 (3") solid
 100 (4") solid with lights
 100 (4") solid
 100 (4") solid with alum. rim
 125 (5") solid
 125 (5") solid with lights
 125 (5") solid with alum. rim
 150 (6") pneumatic
 150 (6") solid


Dealer

Push handles (17)


Standard




Height adjustable




Abducted height adjustable




Curved rear frame
 (it is not possible to assemble the small wheels for narrow passages)



With straight rear frame (19)



With small wheels for narrow passages



(need technical approval)

Unbalancing system (20)




Right
 Left



Right
 Left


(Only with curved rear frame)

Anti-tip device



Right
 Left
 Right & left

Revolving anti-tip device
 (Only with curved rear frame)




Right
 Left
 Right & left

Single drive (21)
 (camber can be only 0°)
 (rear wheels can be only 22" or 24")




Right
 Left

Extended wheel plate (22)



Head-rest (23)

Padded



Height adjustable

Orbital


Anatomic



Moulded




With lateral supports



Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = cm d = cm


Abductor (24)



1 3 5 7 cm
 cm

Choose or specify the height from the seat cover

Adjustable lateral supports (25)



Right
 Left

Specify the height "h" from seat cover to the center of the swing away


Anodized set (27)
 (Rear wheel rims & hub if Light and forks)

- Black
- Silver
- Light blue
- Red
- Orange
- Cyclamen light purple
 (Light wheels only)

Table (26)

- Trasparent
 - Soft padded
- Specify the fastening:
 single double
 (Only with arm-rest)

Lined cushion (28)
 in foam polyurethane



- 3 cm high
- 5 cm high

Brakes (29)

- Classic - short handle (from P32cm)
- Classic - pull to brake
- Classic long handle
- Brakes inside the sideguard*
- Right single brake
- Left single brake
- Drum
- Drum + classic

* compulsory for some configurations

Options

- 30 "M" velcro (hook) on seat
- 31 "M" velcro for VICAIR cushion
- 32 "F" (loop) on seat
- 33 Seat/backrest connection flap
- 34 Padded rigid cushion
- 35 Chest harness
- 36 Chest belt
- 37 Abductor belt
- 38 Double cross bar *
- 39 Stretched bar
- 40 Pump for pneumatics
- 41 45° pelvic belt
- 42 45° padded pelvic belt
- 43 Black backpack
- 44 Front protective covers

* need technical approval

Notes