

Dealer _____

Key standard option

CHILDREN 3000
 Fold.alum.frame for Kid


CHILDREN 3000 24* 26* 28* 30 32 34 36

L Seat width * Not available with detachable foot-rests

CHILDREN 3000 PLUS

L Seat width (adjustable within the selected range) RANGE 1 24* 26* 28* 30 32 34
 RANGE 2 26* 28* 30 32 34 36 38
* Not available with detachable foot-rests

P Seat depth 28 30 (adjustable from 28 to 34 cm)
32 34 36 38 (adjustable from 32 to 38 cm)
 HS Back-rest height 25 27 29 31 33 35 37 39 41
 HA Front seat/floor height 40 42 44 46 48
 HP Rear seat/floor height 36 38 40 42 44 46
 DP Foot-rest/seat distance 17* 19* 21* 23 25 27 29 31 33
35 37 39
* Not available with detachable foot-rests (unless the front frame is abducted)



cm

Frame colour ①

- Deep sky light blue Orange Pearl red Metallic pink
- Neon orange Neon green
- (see "Colour" catalogue)

Backrest upholstery ②

- Tension adjustable nylon back-rest
- Tension adj. breathable b.
- Without back-rest
- Without seat

Upholstery colour (only nylon) ③

- Black Light blue
- Red Green
- Cornflower blue Orange

Front frame ④

- Standard (Not suggested for L < 32cm)
- Abducted

Foot-rest ⑤

- Detachable
- Not detachable
- Elevating (detachable) Subject to technical approval

Separate aluminium profile footplates ⑥

-
- Foot-holding belts

Single aluminium profile footplate ⑦

-
- Foot-holding belts


Side-guards ⑧

-
- Plastic (standard)** (Not available with reclining backrest)
- Aluminium (curved, same colour)


Height adjustable "T" armrest ⑨

-
- Elbow-rest**
- Classic (23 cm long)
- Classic (30 cm long)
- Straight (26 cm long)


Standard push handles



Straight height adjustable push handles



Abducted height adjustable push handles ⑩



Rear wheels ⑪

- 24" Light (24 spokes) (*)
- 22" Light (24 spokes)
- 20" Light (20 spokes)
- Tyres**
- Pneum. high pressure smooth
- Pneum. high pressure plugged
- Solid Grey 20"/22"x1/24x1.3/8
- Solid black KIK 24" x1
- Marathon Plus 22" or 24"
- (*) Only with P = 34 or 36 cm and HP from 40 to 46 cm

Handrims ⑫

- Aluminium silver
- Chromium-plated steel
- Dark aluminium
- Rubber coated
- Carbolife CURVE L 24"
- Carbolife GEKKO 24"
- Silicon covers

Spoke guards ⑬

- Transparent
- Black
- Fantasy
- code pattern:
- (see catalogue of spokeguards)

Front castors ⑭

- 80 (3") solid with lights
- 80 (3") solid
- 100 (4") solid with lights
- 100 (4") solid
- 100 (4") solid with alum. rim
- 125 (5") solid
- 125 (5") solid with lights
- 125 (5") solid with alum. rim
- 150 (6") pneumatic
- 150 (6") solid

Rear wheel camber

- 1° 3° 5°

- Quick release axles Ø12mm with 42mm bushes (def.)
- Quick release axles Ø1/2" (12,7mm) with 42mm bushes
- Quick release axles Ø1/2" (12,7mm) with 50mm bushes

⑮

