


Dealer

Key standard option

QUASAR ADJ FOR DANCING

QUASAR ADJ

- L** Seat width 34 36 38 40 42 44
P Seat depth 36 38 40 42 44 46
HS Back-rest height 27 29 31 33 35 37 39 41
HA Front seat/floor height 46 48 50
HP Rear seat/floor height 40 42 44 46 48
DP Foot-rest/seat distance 37 39 41 43 45 47 49 cm



Frame colour ①

- Shiny titanium Glossy black
 Pearl white Matte black
 Natural titanium
 (see "Colours" catalogue)
 Neon green Neon orange

Upholstery ②

- Back-rest:**
 Tension adj nylon
 Breathable tension adj
 Without back-rest:
 FAMEKH
 Other seating system.....
Seat:
 Fix (standard, screws)
 Tension adj. (velcro)
 Carbon fibre

Upholstery colour ③

- Black
 Coloured (nylon).....
colours: orange, red, green, cornflower, light blue
 Coloured (breathable).....
colours: red, dry green, grey


Single tubular titanium footplate ④

- Open Carbon fibre closed Tilt adjustable Carbon fibre closed
-

Push handles ⑤


- Without pushing handles Fixed pushing handles Foldable pushing handles
-

Distance between the hub centre and back-rest ⑥




- 9 cm (standard) 7 cm
 0 cm (biamputated) 11 cm
 3 cm 13 cm
 5 cm cm

Backrest tilt ⑦



- 89°
 91° (standard)
 93°

Front frame angle ⑧



- 83° (standard) (A = 16 cm)
 75° (A = 21 cm)

Rear wheel camber ⑨


- 1° 2°
 3° 4°
cm (max 12°)

- Quick rel. axles Ø 12 mm (default)
 Quick rel. axles Ø 1/2" (12,7 mm)
 Cambered bushes length 42 mm (default)
 Cambered bushes length 50 mm

Side-guards ⑩

- Aluminium with carbon fibre wing Detachable
 Aluminium with straight carbon fibre Detachable
-

Backrest bar ⑪



- Standard: P(bar) 7,5cm; H(bar) 19cm; dia 19mm.
 Custom: P.....cm; H.....cm;

Brakes ⑫

- Without any brake
 Scissor brake
 Composite
 Standard
 Pull brake
-

Central revolving anti-tip device (stand.) ⑬

- Small wheel of anti-tip with lights
 Reinforced central antitip
 Reinforced central antitip with carbon footplate for 2nd dancer

Notes

Dealer _____

Rear wheels (14)

- 24" Ultralight (24 spokes)
- 24" Light (24 spokes)
- 24" Star 6x3 (18 spokes, silver or black)
- 24" Spinerger
spokes colour.....(black stand.)
- 25" Spinerger
spokes colour.....(black stand.)
- 25" Light (24 spokes) black

Tyres

- Pneum.Schwalbe One 24"/25"
- Pneum. high pressure Smooth
- Pn.Schwalbe Right Run 24"
- Solid KIK 24" not available for Ultralight wheels

Pushrims (15)

- Aluminium silver
- Aluminium dark
- Chromium-plated steel 24"
- Rubber coated
- Titanium
- CarboLife CURVE L 24" *
- CL Curve L TetraGrip 24" *
- CL Quadro rubber ctd 24" *
- CarboLife GEKKO 24" *
- Silicon covers 24"

* CL pushrims: only GEKKO suitable for Star 6x3 wheels

Spokes guards (16)

- Transparent 24"
 - Black 24"
 - Fantasy
code pattern
- (see catalogue of spokeguards)

Lined cushion
 in foam polyurethan (18)




- 3 cm high
- 5 cm high


Front castors (17)

- 80 (3") solid
- 80 (3") solid with lights
- 100 (4"x1.1/4) solid with plastic rim
- 100 (4"x1) solid with lights
- 100 (4"x1.1/4) solid with alum. rim
- 125 (5"x1.1/4) solid with plastic rim
- 125 (5"x1.1/4) solid with alum. rim
- 125 (5"x1) solid with lights
- 150 (6") pneumatic
- 150 (6") solid

Pockets (19)



- Small under the fix seat
(for fix seat by screws)



- On the back-rest

Anodized set (20)

(Rear wheel rims & hubs if LIGHT 24", sideguard supports, forks and fittings)


- Black
- Silver
- Light blue
- Red
- Orange
- Cyclamen light purple

(only 24" Light wheels)


Options

- 22 "M+F" velcro strips (hook and loop) on seat
- 23 "M" velcro strips (hook) on seat
- 24 "F" velcro strips (loop) on seat
- 25 Rigid welded seat
- 26 Welded sideguards, in titanium
- 27 Pump for pneumatics
- 28 45° pelvic belt
- 29 45° padded pelvic belt
- 30 Black backpack
- 31 Front protective covers


STANDARD SEAT




ERGO SEAT




"Y" SHAPED FRONT FRAME
 (stand.)



"V" SHAPED FRONT FRAME



STRAIGHT FRONT FRAME



Custom front opening
cm

(The tilt adjustable footplate only is suitable with the "V" shaped front frame)

(Standard front frame for seat width of 34cm)