

Dealer



Key standard option

- ①
- STAND-UP ST** **STAND-UP ST JUNIOR**
 Static standing with powered elevation for upright position
- STAND-UP EE** **STAND-UP EE JUNIOR**
 Manual propelled standing with powered elevation for upright position
- STAND-UP ETM**
 Powered elevation and traction standing for upright position

Seat width ②

- 36 cm (Stand-Up EE, ST Junior)
 40 cm (standard)
 44 cm
 46 cm



Thoracic support height ③

- 24 cm 
 29 cm 

(*) not available on Stand-up Junior

User measurements ④

A - User height cm
 B - Beneath kneecap height cm



(B: standard 45 cm)

Frame colour ⑤

- Metallic grey
 Metallic blue
 Deep sky light blue

 (see "Colours" catalogue)

Options ⑥

- Pelvis supports

STAND-UP

STAND-UP ST
"with powered elevation"



STAND-UP EE
"with powered elevation and manual traction"



STAND-UP ETM
"with powered elevation and traction"

