

Order N° \_\_\_\_\_  Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer Ref. \_\_\_\_\_ User's weight Kg \_\_\_\_\_

PAG. 1/2  
Jan 2022




Dealer

**Key**     standard     option

# ALHENA

- L** Seat width    34 36 38 40 42 44 46  
**P** Seat depth    34 36 38 40 42 44 46  
**HS** Back-rest height    29 31 33 35 37 39 41  
**HA** Front seat/floor height    42 43 44 46 48 50  
**HP** Rear seat/floor height    38 40 42 44 46  
**DP** Foot-rest/seat distance    33 35 37 39 41 43 45 47    cm



- Frame colour**  
 Shiny titanium     Matte black  
 Pearl white     Natural titanium  
 Neon green     Neon orange  
 .....  
 (see "Colours" catalogue)

- Upholstery**  
 Tension adjustable nylon back-rest  
 Tens.adjustable breathable backrest  
 Without backrest  
 **FAMEKH** custom postural backrest  
(fill up the proper order form)

- Upholstery colour**  
 Black  
 Coloured.....  
 Available colours:  
 NYLON: orange, red, green, cornflower blue, light blue;  
 BREATHABLE: dry green, red, grey;

- Not detachable foot-rest** ④  
 (with DP=43 cm and P=40 cm)  
 (standard)  
  
 84°  
 16 cm

- Front Frame** ⑤  
 STANDARD  
 Y  
(only in case of not detachable foot-rest)  
 V  
(only in case of not detachable foot-rest & without any automatic closure footplate)

- Closed carbon single foot-plate** ⑥  
  
 Up & outward  
 Foot-holding belts

- Automatic closure foot-plate in carbon** ⑦  
  
 Foot-holding belts

- Not detachable foot-rest**  
 (with DP=43 cm and P=40 cm)  
  
 76°  
 21 cm  
 **Detachable foot-rest**  
 84° (std)  
 76°  
 (approx.measurements)

- Separate aluminium profile tilt adjustable foot-plates** ⑧  
 forward     backward  
  
 Foot-holding belts

- Single aluminium profile tilt adjustable foot-plate** ⑨  
 forward     backward  
  
 up & outward  
 Foot-holding belts

- Titanium tube structure** ⑩  
  
 open  
 closed by carbon plate  
 Forward     Backward

- Side guards** ⑩  
 Aluminium with carbon fibre wing  
 Plastic  
 Curved aluminium (same colour or black)  
 Straight aluminium (same colour or black)

- Tip-up, detachable sport arm-rest h. 22 cm** ⑪  
  
**Elbow-rest**  
 Tubular  
 Classic (23 cm long)

- Tip-up height adjustable from 22 to 28 cm "L" arm-rest** ⑫  
(detachable, open with sideguards)  
 Not detachable  
  
**Elbow-rest**  
 Tubular  
 Straight (26 cm long)  
 Classic (23 cm long)  
 Classic (30 cm long)

- Rear wheels** ⑬  
 24" Ultralight (24 spokes) black  
 24" Light (24 spokes)  
 24" Star 6x3 (18 sps) silver or black  
 24" (36 spokes)  
 24" Mountain-bike CITYBIKE tyre  
 24" Mountain-bike OFFROAD tyre  
 24" Spinerger  
 25" Spinerger  
 25" Light (24 spokes) black  
**Tyres**  
 Pneumatic Marathon Plus  
 Pneum. high press.plugged 24"  
 Pneum. high press. smooth  
 Solid 24"(not avail.with Ultralight)

- Pushrims** ⑭  
 Aluminium silver  
 Aluminium dark  
 Chromium-plated steel 24"  
 Rubber coated  
 Titanium  
 Carbolife CURVE L24"\*  
 CL Curve L TetraGrip 24"\*  
 CL Quadro rubber ctd 24"\*  
 Carbolife GEKKO 24"\*  
 Silicon covers  
 \* not suitable for wheels Star 6x3 and Mountain Bike