


Order N° \_\_\_\_\_  Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer Ref. \_\_\_\_\_ User's weight Kg \_\_\_\_\_

**PAGE**  
Jan 2022



Dealer

**Key**       standard       option



- L** Seat width      34 36 38 40 42 44 46
- P** Seat depth      36 38 40 42 44
- HS** Back-rest height 27 29 31 33 35 37 39
- HA** Front seat/floor height 46 48 50
- HP** Rear seat/floor height 39 41 43 45 47
- DP** Foot-rest/seat distance 31 33 35 37 39 41 43 45 47 49 cm

- Frame colour**
- Shiny titanium
  - Pearl white
  - Natural titanium
  - ..... (see "Colours" catalogue)
  - Neon green
  - Glossy black
  - Matte black
  - Neon orange

- Upholstery**
- Tension adjustable nylon back-rest
  - Tension adj. breathable back-rest
  - Without back-rest

- Upholstery colour**
- Black
  - Coloured (nylon).....  
colours: orange, red, green, cornflower, light blue
  - Coloured (breathable).....  
colours: red, dry green, grey)

**Single tubular foot-plate**

- Open
- Carbon fibre closed
- Aluminium closed
- Tilt adjustable Carbon fibre closed
- Footrest colour: black (standard)
- Foot-rest colour.....

**Pushing handles**

- Without pushing handles
- Standard pushing handles
- Height adjustable pushing handles
- Foldable pushing handles

**Distance between the hub centre and back-rest**

- 4 cm
- 5 cm
- 7 cm (standard)
- 9 cm
- 11 cm

**Back-rest tilt**

- 89°
- 91° (standard)
- 93°

**Seat with shock absorbers**

- User weight : ..... kg
- It is not possible to have two anti-tip devices

**Side guards**

- Aluminium with fixed carbon fibre wing
- Straight aluminium
- Aluminium with detachable carbon fibre wing

**Detachable "L" arm-rest - h 22 cm**

- Tip-up Elbow-rest
- Tubular
- Straight (26 cm long)
- Classic (23 cm long)
- Classic (30 cm long)

**Brakes**

- Push brake (stand.)
- Pull brake
- Composite brake
- Scissor brake

**Revolving anti-tip device**

- Right
- Left

Notes

