

QPX

**TELAIO IBRIDO
IN TITANIO
O ALLUMINIO**

TITANIUM OR
ALUMINIUM
HYBRID FRAME



AGILE COME UN TELAIO RIGIDO, PRATICA COME UNA PIEGHEVOLE L'innovativa crociera, da cui prende vita il modello, è progettata e realizzata con la massima precisione. Garantisce stabilità e sostegno, scomparendo alla vista.

PERFETTA PER OGNI BAGAGLIAIO

L'azione di apertura/chiusura impiega pochi secondi e minimizza l'ingombro.

IL TEMPO È PREZIOSO QPX lo compatta.


THE AGILITY OF A RIGID FRAME, THE VERSATILITY OF A FOLDABLE WHEELCHAIR The innovative crossbar design, which gives life to the wheelchair is machined with clockwork precision. It grants sustain and stability, disappearing from view.








FITS EVERY TRUNK The opening and closing actions require mere seconds and the space taken is minimal.

TIME IS VALUABLE have some more, with QPX.

Rigido +
Pieghevole =
Ibrido

Rigid +
Folding =
Hybrid



								
	1° 4°				27			
34	52 58		A: ~ 41	36	29			seduta seat
36	54 60			38	31	HA: 46÷50		38x40 cm
38	56 62	25		40	33		P = 40 cm α = 83°	
40	58 64		B: ~ 80	42	35	HP: 40÷48	~ 84 cm	10,3 kg (Ti)
42	60 66			44	37			10,6 kg (Al)
44	62 68			46	39			
					41			

Portata massima 120 kg / Max user weight 120 kg

Alcune misure sono approssimate in base alle medie / Some measurements are averaged



Crociera precisa e
innovativa
Precise revolutionary
crossbar



Salva abiti fissi o
rimovibili in carbonio
Fixed or removable carbon
clothes guards



Baricentro regolabile
Adjustable centre of
gravity



Pedana automatica in
carbonio
Self closing carbon
footrest



Minimo ingombro
(chiusa)
Minimum width