

ALHENA

**TELAIO
PIEGHEVOLE
IN TITANIO**

TITANIUM
FOLDING
FRAME



ELEGANZA E LEGGIADRIA NEL DNA Carrozzina superleggera pieghevole, interamente in titanio. Unica e raffinata, coniuga forma e piacevolezza d'uso.


COMPLESSA SEMPLICITÀ ALHENA è leggera, sofisticata e versatile. La piastra HWP (Helicoid wheel Plate), oltre a proporre un'ampia scelta di posizioni per la ruota posteriore, dona carattere al design del telaio.

FORM AND FUNCTION Ultralight foldable wheelchair, entirely in titanium. Unique and refined, it combines beauty and ease-of-use.

SIMPLY COMPLEX ALHENA is lightweight and versatile, while its presence is elegant, its design polished. The HWP plate (Helicoid wheel Plate) brings the whole design together, as well as being a tremendous setup tool.

Fiera leggerezza

Ultralight prowess



	1° 4°		34	29			
34	52 57		36	31		P = 40 cm COG = 6 cm	seduta seat
36	54 59		38	33	HA: 42÷50	pedana fissa fixed footrest ~ 87 cm	38x40 cm
38	56 61	29	40	35	HP: 38÷46	p. estraibile detachable foot. ~ 95 cm	9,7 kg
40	58 64		42	37			
42	60 66		44	39			
44	62 67		46	41			

Portata massima 120 kg / Max user weight 120 kg

Alcune misure sono approssimate in base alle medie / Some measurements are averaged



Pedane estraibili (opt.)
Removable footrests (opt.)



Due differenti forme anteriori del telaio: dritto, a Y (in foto)
Two different front frame shapes: straight, Y shape (pic.)



Pedana automatica in carbonio (opt.)
Self closing carbon footrest (opt.)



Minimo ingombro (piegata)
Minimum volume (folded)