

S O L

TELAIO FISSO
IN TITANIO

TITANIUM
RIGID FRAME



COME UN ABITO SU MISURA Il telaio fisso monotubo, in titanio, è ampiamente personalizzabile in tutte le sue geometrie. Ogni carrozzina è unica, pensata e costruita attorno a una persona specifica.


BELLEZZA ULTRALEGGERA I dettagli in carbonio alleggeriscono e valorizzano la carrozzina, esaltandone agilità e reattività: una vera ultraleggera che fa la differenza ad ogni spinta.

FITS LIKE A GLOVE The open titanium frame is built to spec, with many possible combinations of shapes and sizes. Each wheelchair is unique from start to finish, built with a single user in mind.

ULTRALIGHT BEAUTY The carbon details decorate the wheelchair and help keeping it light and responsive: you can feel the magic at every push.

Essenziale Originale Speciale

Essential
Outstanding
Steadfast



	cm	L	cm	cm	A	B	P	cm	HS	cm	HP	HA	cm	a	COG	asse fisso	asse regolabile
		1°	4°														
34	54	60		43			36			27						seduta	seat
36	56	62		45			38			29						38x40 cm	
38	58	64		47			40			31						7,5 kg	
40	60	66		49			42			33						asse fisso	fixed COG
42	62	68		51			44			35						~ 82 cm	
44	64	70		55			46			37						7,8 kg	
46	66	72		57						39						asse regolabile	adjustable COG

Portata massima 120 kg / Max user weight 120 kg

Alcune misure sono approssimate in base alle medie / Some measurements are averaged



Baricentro fisso o regolabile
Fixed or adjustable centre of gravity.

Freni a scomparsa
Scissor brakes

Sedile ergonomico o standard
Ergonomic or standard seat

Pedana in carbonio regolabile
in inclinazione (opt.)
Tilt adjustable carbon footplate (opt.)