

QUASAR

**TELAIO FISSO
IN TITANIO**

**TITANIUM
RIGID FRAME**

ULTRALEGGERA PER ANTONOMASIA

Ogni elemento è studiato per minimizzare il peso, con una configurazione completa (spondine, ruote e freni) a partire da 7 kg.

INEDITA QUASAR è completamente personalizzabile. Ogni singolo telaio è disegnato e realizzato per fondersi con l'utente valorizzandolo. Dinamica, è un vero passepartout per tutte le occasioni della vita.

QUINTESSENTIALLY ULTRALIGHT

Every component is designed to minimize weight, with a complete setup (wheels, sideguards and brakes) starting at 7 kg.

UNPRECEDENTED QUASAR is entirely bespoke. Each frame is designed and built to enhance the user, seamlessly, and to make you feel at ease in every moment, up to any challenge.



Essenzialmente ultraleggera

Ultralight to the core



cm L	cm	cm	A B cm	P cm	cm HS	HP cm	HA cm	cm
	1°	4°						
34	52	58	44					seduta seat
36	54	60	46	A: ~40	36	27		38x40 cm
38	56	62	48		38	29		7 kg asse fisso fixed COG
40	58	64	50		40	31	HA: 46÷50	
42	60	66	52	B: ~78	42	33		P = 40 cm α = 83° COG = 9 cm
44	62	68	54		44	35	HP: 40÷48	~ 78 cm
46	64	70	56		46	37		7,4 kg asse regolabile adjustable COG
						41		

Portata massima 120 kg / Max user weight 120 kg

Alcune misure sono approssimate in base alle medie / Some measurements are averaged



Baricentro fisso o regolabile
Adjustable or fixed centre of gravity.



Tre forme anteriori differenti:
dritto, a V, a Y (in foto)
Three front frame shapes:
straight, V shape, Y shape (pic.)



Sedile ergonomico o standard
Ergonomic or standard seat



Minimo ingombro
Minimum volume